



Orien
Speak with your eyes



Scan with the iPhone
camera to get Orien.

Getting started with Orien

A free iPhone app that lets someone talk using only their eyes. A setup guide for families and caregivers.

What you need

- An iPhone with Face ID and iOS 17 or later (iPhone XR, XS or newer). iPhone SE won't work, and iPad isn't supported yet.
- A mount or stand that holds the phone steady.
- Good light on the face.
- A charger nearby. The screen stays on while Orien is open.

How to talk with Orien

- Look up, down, left or right to highlight a choice.
- Hold your gaze until the bar fills, or blink on purpose with both eyes.
- Look down to go back. On the first page, down opens People.
- Look right at Next Page for more choices.

Set it up

- ① **Install Orien.** Scan the code above, or search Orien Speak with Eyes in the App Store. It's free, with no account.
- ② **Mount the phone** in front of the face, at eye level, about as far away as for a FaceTime call.
- ③ **Open Orien** and tap Allow for the camera. Contacts and Photos are optional.
- ④ **Follow the short tour.** When OK appears, look at it or tap it, then enter a name. Tapping the keys is fine.
- ⑤ **Try it.** Look up to open E-Talk, then choose "I'm in pain."

Too slow or too fast? Look right twice to reach the third page, look left for Settings, then choose Speed.

For caregivers

- You can tap anything on the screen. A tapped phrase is spoken right away.
- If "Face not detected" appears, center the face in front of the phone, add light, and make sure nothing covers the top of the phone.
- No sound? Turn the volume up and make sure silent mode is off.
- Texts open in Messages and need a tap to send.

Orien is not for emergencies. In an emergency, call 911.

Given to you by _____

Phone _____

Help and questions

speakwithorien.com/support
speakwithorien.com/contact